

FXE Instrument Parameters for Early User Operation

All parameters are subject to change, pending the commissioning process of both the accelerator and the instrument.

Please discuss your experiment plans with your FXE instrument scientist **before** submitting your proposal. They can help you with any details that may have been updated, assist with evaluating the experiment feasibility, and much more.

Contacts: fxe-instrument@xfel.eu

sample.environment@xfel.eu

useroffice@xfel.eu

X-Ray Photon Beam Parameters

| | | |
|------------------------------------|--|---|
| photon energy | 9 keV | Possible option: adjustable 7 -14 keV |
| pulse energy | 0.5 mJ | at 9 keV |
| photons per pulse | 3.5×10^{11} | at 9 keV |
| pulse duration | 50 fs | calculated from e- beam properties |
| spot size on sample | ~ 10 μm (focus) | variable up to ~100 μm |
| photons/ μm^2 on sample | $>3.5 \times 10^9$ | derived |
| train repetition rate | 10 Hz | fixed |
| intra-train repetition rate | 1.1 MHz | Possible option: 4.5 MHz, possibly 100 kHz, arbitrary pulse pattern |
| no. of bunches per train | 1-60 | Possible option: from 1 to 300 |
| pointing stability | 2 μrad | possibly drifting over entrance apertures |
| $\Delta E/E$ | ~ 0.2% | calculated |

Scattering: Large Pixel Detector Parameters

| | | |
|--------------------------|---|---|
| no. of pixels | 1024 × 1024 | 4 quadrants, each 512 × 512 pixels |
| pixel size | 500 × 500 μm^2 | |
| sensor | Si, 500 μm | |
| distance sample-detector | 100 mm – 1.5 m | on motorized stage, up to 8 m manually |
| central tube diameter | 4, 10 mm | Q_{min} (9 keV, 100 mm distance) = 0.13 \AA^{-1} , 0.32 \AA^{-1} |
| central hole size | 1 – 30 mm | |
| max. Q range at 9 keV | ~ 5.1 \AA^{-1} | at the edge of detector for hole-centred beam (100 mm) |
| dynamic range | 10^5 at 12 keV | |
| quantum efficiency | 89% at 12 keV | 98% at 10 keV, 38% at 20 keV |

Possible Added Instrumentation*

| | | |
|---------------------------------------|---------------------------|---------------------------------------|
| Laser-X-ray timing tool | < 100 fs | Cross correlation time |
| Single shot spectrum analyzer | ~0.1 eV resolution | at 7 keV, single burst XANES possible |
| 4-bounce primary monochromator | Si(111) | 5-14 keV possible |

*) inquire to FXE staff about the availability


XES: von Hamos Spectrometer

| | |
|---------------------------|--|
| crystal bending radius | 500 mm |
| total number of tiles | 16 (30 × 110 mm ²) |
| angular range (1 crystal) | ~ 2.8° ($\theta = 80^\circ$) ~ 1.15° ($\theta = 60^\circ$) |

Detectors

| | “GreatEyes” (2D) | GOTTHARD (1D) |
|----------------------------|--|---|
| detector efficiency | ~ 90% (5 keV), > 25% (10 keV) | > 90% (5 – 10 keV), > 50% (15 keV) |
| detector frame rate | 10 Hz | 40 kHz – 0.8 MHz |
| no. of pixels (pixel size) | 1024 × 256 (26 × 26 μm^2) | 1280 × 1 (50 μm × 8 mm) |
| pixel dynamic range | ~ 350 at 7 keV | 10^4 at 12 keV |

Available Crystals and Energy Coverage in 60° – 90° Geometry



XES: Johann Spectrometer

| | |
|------------------------------|---|
| crystal radius of curvature | 1 m |
| total number of crystals | up to 5 10 cm diameter |
| detector and characteristics | APD, PIN diode MHz readout possible |
| | GOTTHARD see above |
| available crystals | Si(220), Si(111), Si(531), Si(400), Si(620), Si(551) (* inquire with FXE staff about the availability) |

Sample Delivery Systems

| | | |
|--------------------------|---|--|
| liquid flat-sheet jet | 100 μm, 300 μm | fixed, for wide range of viscosities |
| liquid flat-sheet jet | 1 – 100 μm | variable, for water and low viscosity liquids (inquire with FXE staff in case of interest) |
| slow fixed target sample | “XFEL design” | |

Optical Laser Systems Parameters

Three **synchronised femtosecond to picosecond laser systems** will be available. In addition, CEP-stable **single-cycle terahertz pulses** (inquire to FXE staff about the availability) can be generated using optical rectification, with a centre frequency of 0.3 THz (1 mm wavelength) and ~50 μJ pulse energy at 100 kHz repetition rate. All laser pulses can be time delayed with respect to the X-ray pulse over a range of 4.6 ns in steps of 2.5 fs. Optical parametric amplification and white light generation schemes are being investigated and will be available in Run 2.

| | pump-probe laser system I | | pump-probe laser system II | | pump-probe laser system “Tangerine” | |
|--|--|--|-----------------------------------|--|---|---|
| wavelength | 800 nm | | 1030 nm | | 1030 nm | |
| pulse duration | 15 – 300 fs | Close to TL | 0.9 or 500 ps | Compressed or chirped | 300 fs | |
| train repetition rate | 10 Hz | burst duration up to 300 μs | 10 Hz | burst duration up to 300 μs | 0.188 – 4.5 MHz | variable, burst mode operation possible |
| intra-train repetition rate | 4.5 MHz or 1.1 MHz | variable, down to single pulse | 4.5 MHz or 1.1 MHz | variable, down to single pulse | 4.5 MHz | variable, down to 100 kHz |
| wavelength conversion | SHG, THG | no OPA in this Run | SHG, THG, THz | | SHG, THG, FHG | no OPA in this Run |
| pulse energy (fundamental) | 50 μJ / 200μJ | at 4.5 / 1.1 MHz | 1 mJ / 4mJ | at 4.5 / 1.1 MHz | ~ 6 μJ / 25 μJ | at 4.5 / 1.1 MHz |
| efficiency of 2 nd harmonic gen | t.b.d. | | t.b.d. | | ~ 40% at 400 kHz | |
| efficiency of 3 rd harmonic gen | t.b.d. | | t.b.d. | | ~ 10% at 400 kHz | |
| efficiency of 4 th harmonic gen | N/A | | t.b.d. | | ~ 5% at 400 kHz | |
| arrival time jitter w.r.t. X-rays | <100 fs | estimated | <100fs | estimated | t.b.d. | |