

Ideas for microfluidics experiments at MID

Sarah Köster, Tim Salditt

Institut für Röntgenphysik, Georg-August-Universität Göttingen

Early Science Workshop @ MID 27.01.2015



Cellular imaging – light microscopy





- live cell imaging
- molecule specific labeling
- extensive sample preparation
- small penetration power
- spatial resolution: ~200 nm (but: STED, STORM!!)

images: Rosmarie Sütterlin, Ueli Aebi (Biozentrum Basel)

Cellular imaging – electron microscopy



SEM image by P. Walther, M. Beil, Uni Ulm



TEM images by W. Möbius, MPI EM Göttingen

- static samples (dry, cryo...)
- imaging of surfaces/slices
- very extensive sample preparation
- spatial resolution: nm

Sarah Köster, University of Göttingen

Cellular imaging – X-rays?!



- Introduction
- Imaging cells with X-rays: nano-diffraction, ptychography
- Sample environments: microfluidics
- Experiments @ FLASH and SACLA & Ideas for microfluidics experiments @ MID

Networks in epithelial cells



keratin bundles



SEM image by P. Walther, M. Beil, Uni Ulm

- occur mainly in epithelial cells
- provide tensile strength
- single filament thickness ≈ 10 nm

Sarah Köster, University of Göttingen

Nano-diffraction setup





- small beam size on sample $\approx 150 \times 150 \text{ nm}^2$
- scattering: high resolution in real space and in reciprocal space
- *local, internal* bundle structure and orientation
- ID 13/ESRF, P10/PETRAIII, cSAXS/SLS

Weinhausen et al, New. J. Phys. 2012, Weinhausen & Köster Lab Chip 2013

Strategy



(1) chemically fixed & plunge frozen & freeze-dried cells

- good electron density contrast
- can be prepared beforehand
- no biosafety-level necessary
- 2 chemically fixed & hydrated (buffer) cells
 - fewer preparation steps
 - "test": aqueous environment
 - preparation beforehand
 - no biosafety-level necessary
- 3 living cells (medium)
 - dynamic measurements (microfluidics)
 - "in situ", "in operando"

X-ray dark field images

fixed, freeze-dried cells

ID13/ESRF

visible light fluorescence (inverted gray scale)





step size 2 μm

Internal nanostructure



fixed, freeze-dried cells



Weinhausen et al, New. J. Phys. 2012

Sarah Köster, University of Göttingen

Hair cell stereocilia

actin bundles

- sensors in the inner ear
- highly parallel bundles of actin
- arrangement of the filaments?



David Furness, Keele University







а

Internal nanostructure



actin; fixed, freeze-dried cells





Scattering geometry









Piazza et al., ACSNano, 2014



Simulation

Tim Salditt, Göttingen







Priebe et al., Biophys. J., 2015

Sarah Köster, University of Göttingen

Ptychograpy





Piazza et al., ACSNano, 2014

Sarah Köster, University of Göttingen



- Introduction
- Imaging cells with X-rays: nano-diffraction, ptychography
- Sample environments: microfluidics
- Experiments @ FLASH and SACLA & Ideas for microfluidics experiments @ MID

A step forward: hydrated cells



measurements on hydrated cells

- native protein structure
- no artifacts from freezing/drying
- lower electron density contrast
- absorption/scattering due to water layer
- increased radical production and mobility in water

sample environment

- compatible with x-ray diffraction
- compatible with cell culture
- small/light enough to fit into the set-up

Sarah Köster, University of Göttingen

A microfluidic chamber for cells



- fluid flow: nutrient supply and waste removal
- Si₃N₄ membrane windows as substrate and window material for X-rays



A microfluidic chamber for cells









• channel size $\approx 500 \times 150 \ \mu m^2$

Fixed-hydrated cells

fixed, hydrated cells



x 10⁵



good signal-to-noise ratio in dark-field image

Radial intensity





- different scan regions
- averaged scattering patterns
- azimuthal integration
- background subtraction
- power-law fit



Initially living cells

living cells

- P10/PETRA III
- transport in microtubes filled with mediumstorage in petri dishes incubator



asymmetric scans to skip damaged regions

Radial intensity









 different (higher) power law exponents for living cells compared to fixed-hydrated cells

Weinhausen et al, PRL 2014

Sarah Köster, University of Göttingen

Structural changes





Conclusions

- X-rays provide a complementary approach for imaging cells (light microscopy, electron microscopy)
- different (contrast) methods are currently being developed and improved
- radiation damage remains a great challenge
- reciprocal space provides structural information at high resolution



Howells et al., 2009

A

- Introduction
- Imaging cells with X-rays: nano-diffraction, ptychography
- Sample environments: microfluidics
- Experiments @ FLASH and SACLA & Ideas for microfluidics experiments @ MID

Sarah Köster, University of Göttingen

Requirements @ MID

- scintillation-based detector (AGIPD) @ ~ 10 m from focus
- possibly second detector, "drilled-through"
- optical microscopy for sample alignment
- feed-through for connections to microfluidic chambers



Acknowledgements

Göttingen:

Tim Salditt Susanne Bauch Britta Weinhausen Oliva Saldanha Clement Hemmonot Valeria Piazza

Wiebke Möbius Tobias Moser







Juliane Reinhardt

Dresden: Jens Patommel

Winheim



Heidelberg: Harald Herrmann, Norbert Mücke, Stefan

Grenoble: Manfred Burghammer, Michael Reynolds

Andreas Menzel

Aachen: Rudolf Leube, Reinhard Windoffer, Anne Kölsch

Hamburg: Anders Madsen, Jörg Hallmann, Christian Schroer,





Excellence Initiative SFB 755, SFB 937, CMNPB, KO 3572/5-1

> Deutsche Forschungsgemeinschaft

Deutsche Telekom Stiftung









Bundesministerium für Bildung und Forschung

DFG Sarah Köster, University of Göttingen