# 11<sup>th</sup> call @ Femtosecond X-ray Experiments

Peter Zalden on behalf of FXE

European XFEL FXE – Femtosecond X-ray Experiments

Last town hall meeting: 12.5.2022

#### Agenda 20.4.2023

- 1. Overview FXE
- 2. News since last town hall
- 3. Questions/discussion

European XFEL





peter.zalden@xfel.eu

### FXE: Femtosecond (Hard) X-ray Experiments

A suite of **Simultaneous** X-ray Tools combined with flexible laser excitation sources

- Single-shot dispersive resonant and non-resonant XES: von Hamos
- Wide(Small)-angle X-ray Scattering: Large Pixel Detector (LPD) and Jungfrau
- Sample motion stack, goniometer for **single-crystal X-ray diffraction** with Jungfrau detector motion using the robot arm
- Tuneable laser excitation covering 1030, 515 nm, 257 nm (1 ps), 800, 400, 266 nm (15 or 50 fs) and an **OPA** (50 fs, 240 nm to 3  $\mu$ m) with **THz** in development (LiNbO<sub>3</sub>, 0.2-0.3 THz)
- **X-ray absorption spectroscopy** (5-20 keV): **scanning** (Si111 *2-bounce* mono) and **single-shot** (Spectrum analyzer)

**Scanning** resonant and non-resonant **XES** (RXES): Johann

European XFEL



"Scientific instrument Femtosecond X-ray Experiments (FXE): instrumentation and baseline experimental capabilities" A. Galler, et al., J. Synch. Rad., 26, 1432 (2019) "Ultrafast X-ray Photochemistry at European XFEL: Capabilities of the Femtosecond X-ray Experiments (FXE) Instrument" D. Khakhulin, et al., Appl. Sci., 10, 995 (2020)

# **FXE Group Members**

Engineering team







Martin Knoll

Paul Frankenberger

Siti Heder

#### Postdocs



<mark>Han Xu</mark>



Xinchao Huang

Doriana Vinci





Hao Wang Diana Bregenholt Jakobsen **European XFEL** 

#### Leading Scientist



Chris Milne

Sharmistha Paul Dutta

Administrative Support

Maria Peter

#### Scientists



Dmitry Khakhulin



Peter Zalden

Detector Scientist

Hazem Yousef



Frederico Alves Lima



Yifeng Jiang



Yohei Uemura

Soon

Mykola Biednov

Fernando Ardana Lamas



2023

Instrument scientist



Recent alumni

Vandana Tiwari, postdoc

https://www.xfel.eu/facility/instruments/fxe/group members/





## Sample environment: Liquid & Solid state experiments



- Vacuum environment (1e-5 mbar)
- X-ray probe in transmission geometry

Up to 15k samples accessible per filling/evacuation Parallel X-ray emission and scattering compatible Diffraction up to 8.8 Å<sup>-1</sup> at 16.5 keV and  $2\theta_{max}$ =63° courtesy of P. Zalden (peter.zalden@xfel.eu)



- Helium environment
- Open on 3 sides (XES, XAS, WAXS compatible)
- Parallel UV-Vis flow loop to monitor sample
- Jet diameter 25-200 µm
- Bragg angle range 67-83°
- WAXS maximum Q up to 10 Å<sup>-1</sup>

courtesy of F. Lima (frederico.lima@xfel.eu)



courtesy of D. Khakhulin (dmitry.khakhulin@xfel.eu)

#### Scientific Scope of FXE: Measuring ultrafast dynamics with hard X-rays



### Ultrafast (bio)chemical dynamics



8970

8980

11.79

11.8

11.81

Energy in I 11.82

11.83

11.84

11.85

keV

8990

50

100

150

VHS energy in a.u.

200

9000

Au foil, run 59

Energy /eV

9010

November 2022

run163+164

0.02

0.01

0.00

-0.01

-0.02 है

-0.03

-0.04

-0.05

odds







p003495\_run157\_XAS\_I0\_CH2A\_I\_CH2B

evens-odds

3.900

3.875

3.850 GH2A/CH2B 3.825 CH2

3.800

3.775

3.750

laserOn

laserOff

XAS @ FXE: Status update Cu(dmp)2, 20 mM, 200 µm liquid jet 0.030 *k*=6 E0



#### Next steps:

0.7

- Investigate other X-ray energies and edges
- Continue detection optimization

European XFEL

courtesy of Y. Uemura (yohei.uemura@xfel.eu)

#### February 2023 New sample environment: **Goniometer in Helium environment**

- Thin film targets
- Ideal for low E<sub>ph</sub> spectroscopy (Ti, V, Cr K edges)
- tr-EXAFS possible
- Simultaneous diffraction 20 pulses/train (94 kHz)













# SFX-type expts@FXE

Validation of SFX detection with Lysozyme

300						
Ċ	· · ·	ap .	0,0			
	D	ď			•	
200 -			0 0 0			
		о <sup>р</sup>		<b>P</b>		
100	٦	ē.	<b>1</b> 0. <b>B</b> . o	or	. 8	i D
		ter c			<b>.</b>	
		5		2 <b>0</b>	ę o	G
pixels		q - <b>Q</b>	K		•	
	ц. С. С.					
	0	P.	.0	0		
-100 -						
		•	- 10 M	0		
	E	5			d <sup>B</sup>	
-200 -		0				
		E		0000		0
-300 +		0				
30	200 200	100	0 pixels	-100	-200	-30



- smSFX proved itself  $\geq$ feasible and efficient at FXE
- Combination of high  $\geq$ rep.-rate and high photon energy (15 keV) allows for solving small molecule structures with high fidelity after a short acquisition time



### Call 11: Contact us with any questions or for more details



Email: <u>christopher.milne@xfel.eu</u> OR <u>fxe-support@xfel.eu</u>

https://www.xfel.eu/ facility/instruments/fxe/